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
Abb. Therapiestütze

Supplement user manual stabiliser

COMPATIBLE WITH THE FOLLOWING MODELS

Non-motorized bicycles (wheel size 16" – 28") with gear hub.

INTENDED USE

Enables the conversion of an existing two-wheeler to a multi-track vehicle. Acquired limitations such as dizziness and/or increasing unsteadiness when cycling due to age or e.g. artificial joints, etc. can be compensated for.

NOTE

Only compatible with non-motorized vehicles.

AVAILABLE SIZES

Stabiliser 16":

Suitable for all commercially available 16" bicycles.

Stabiliser 20":

Suitable for all commercially available 20" bicycles.

Stabiliser 24":

Suitable for all commercially available 24" bicycles.

Stabiliser 26":

Suitable for all commercially available 26" bicycles.

Stabiliser 28":

Suitable for all commercially available 28" bicycles.

SCOPE OF DELIVERY

Pre-assembled stabiliser

INSTALLATION

1. Remove the therapy support from the packaging and take off the support strut that is fastened with cable ties.
2. Then loosen the adjustment brackets.
3. Mount the adjustment brackets on the axle of the gear hub. Disassemble the adjustment brackets on the stabiliser (A) and then loosen the axle nuts of the gear hub. Place the adjustment brackets on the axle and screw them back on with the required torque.
4. Reconnect the stabiliser and the adjustment brackets and screw them tight. Loosen the sleeve screw at the upper end of the spring element. Insert the spring element into the damper mount and connect both with the sleeve screw.
5. Connect the support strut to the flap on the upper arch of the stabiliser frame. The other end of the support strut is connected to the

luggage carrier mount on the bicycle frame. An adapter may be necessary for this. This adapter depends on the bicycle and is not included in the scope of delivery. If you are unsure about this, we recommend that you consult a specialist.

ADJUSTMENT

1. Loosen the screw connections on the shock mount. Hold the bicycle in a vertical position and adjust the outrigger of the stabiliser so that the wheel touches the ground. Proceed analogously on the other side.
2. Loosen the screw connection of the support strut with the flap on the frame of the stabiliser and bring the stabiliser into a vertical position. Tighten the connecting screw again.

⚠ DANGER

Ensure that the stabiliser does not collide with the riding surface or parts of the bicycle.

⚠ DANGER

Check the screw connections regularly and retighten them, if necessary.

⚠ DANGER

The service life of the accessories depends on many factors, such as the conditions of use or the physique of the rider. Other factors, e.g. falls, improper use can significantly reduce the service life.

⚠ DANGER

Mounting the flexible training wheels is not comparable with mounting training wheels on children's bicycles.

⚠ DANGER

Cross railway tracks and sewer covers carefully to avoid falling. Cross the railway tracks at a right angle if possible. **Danger of accident!**

⚠ DANGER

Let the assembly and adjustment work of the stabiliser be carried out by a specialist. Specialist knowledge and special tools are required for all adjustment, care and maintenance work on the vehicle. **Danger of accident!**

NOTE

Make a test-ride with your new stabiliser.

NOTE

Any modifications to the therapy support are not permitted and will invalidate the warranty.

DIRECTIONS OF USE

Carry out driving exercises with your specialist dealer to familiarise yourself with driving the stabiliser and the new driving experience that comes with it.


Fig. Adjustment of adjustment brackets at the axle nuts of the gear hub

Fig. Stabiliser